



**Purpose:** To obtain an IR scan of pin coatings, bowling ball cover stocks, ball cleaners, lane cleaners, and lane oil.

**Materials:**

- Dell Optiplex GX620
- Spectrum 100 Series Perkin Elmer FT-IR
- Kimwipes
- Tweezers
- Isopropyl alcohol (IPA)

**Procedure:**

1. Power on the Instrument
  - a. Power on FT-IR and wait for display on PE Spectrometer
  - b. Power up computer and log in as Administrator (no password)
  - c. Select the Spectrum icon.
  - d. At Perkin Elmer Log in click OK
  - e. At Spectrum 100 Log in click OK
2. Run a System Suitability Check
  - a. Perform the system suitability check at the beginning of each workday the instrument will be used.
  - b. Using a Kimwipe, swab crystal (located on the top-plate) with IPA and allow the top-plate to dry.
  - c. From the menu select: Instrument > Validate > System Suitability
  - d. Select OK to run system suitability check. Once instrument has been validated, print the report and proceed to Section C.
3. Take a Background Scan
  - a. From the menu select: Instrument > Scan
  - b. In the *Scan and Instrument Setup* window, select the *Sample* tab and enter “background” for the name.
  - c. Select the *Scan* tab and under the options heading, select background from the drop-down menu.
  - d. On the same tab and under the durations heading, enter 15 for the number of scans and click start.
  - e. Begin the background scan by either pressing the Scan button on the screen. Alternately, the blue scan button on the instrument may be pressed to begin the scan.
  - f. Overwrite the previous background scan by saving the newly collected file.
4. Take a Sample Scan
  - a. From the menu select: Instrument > Scan
  - b. In the *Scan and Instrument Setup* window, select the *Sample* tab and log sample number & name under details. Enter all other known information.
  - c. Select the *Scan* tab and under the durations heading, enter 15 for the number of scans.
  - d. Using tweezers place the sample over the crystal and select start.
  - e. Select the appropriate foot, move arm over the sample and apply light pressure by rotating knob clockwise. Monitor the force gauge and do not exceed 150.



- f. Select scan.
  - g. When complete, sample spectrum will appear on monitor.
  - h. Save the collected scan to C:/pel\_data/spectra
  - i. Turn the knob counter clockwise to release the sample and return arm to its original position.
  - j. Remove the sample with the tweezers.
  - k. Clean the top-plate with IPA before analyzing another sample.
5. Compare Spectra (when required)
- a. Load spectral file to compare.
  - b. From the menu select: Setup > Compare > Type
  - c. In the directories heading, select previous file path and choose Delete.
  - d. Select Add and choose the folder containing the data to be compared.
  - e. Select OK.
  - f. From the menu select: Process > Compare OR select the Compare icon located on the toolbar.
  - g. If desired, select the Comparison window and print results.